



EMS Vacuum Immobilisation Catalogue 2023



RedVac EMS Vacuum Immobilisation

Content	Page
RedVac made in Austria	3
The RedVac System	4-5
Mattresses	6-10
Head Immobilisation	11
Extrication Device	11
Extremity Splints	12-18
Pediatric	19-20
Pumps, Bags & Accessories	21
Thermo & HRS Rescue Sacks	22
Helicopter Rescue Sack	23



Developed and manufactured in Austria

Kohlbrat & Bunz products are developed and manufactured 900 meters above sea level in the middle of the Austrian Alps in its own facilities by its experienced and meticulously working staff.

People living in the Alps know very well that products used under the extreme conditions of alpine rescue have to work perfectly.


To meet these rigorous demands our products are developed in close cooperation with the mountain rescue teams and other demanding emergency services.



This cooperation has brought about a selection of the best materials and optimal product design leading to products of outstanding performance and quality. To insure the best possible quality Kohlbrat & Bunz work in accordance with EN ISO 9001 & EN ISO 13485.

Worldwide in use

We are very proud of the confidence shown to us and our products by many well-known and demanding emergency services worldwide and are happy to show some examples of their beautiful logos that they have asked us to put on the RedVac products we have made for them.



The RedVac System

for **PATIENT SAFETY, PERFECT HYGIENE & EASE OF USE**

- **The RedVac CLICK-System**
for attachment of straps and handles

The unique CLICK-System makes removal of straps and handles for repositioning, replacement, cleaning and disinfection very easy!



- **The RedVac Multi-Chamber-Systems**
- provides unmatched **EASE OF USE**

The Multi-Chamber-Systems keeps the internal beads evenly distributed eliminating the need for difficult and time consuming manual spreading of the beads. This enables fast, safe and easy application in all situations.

The RedVac Multi-Chamber-System!

The internal sack system keeps the beads evenly distributed in up to 120 chambers.



The RedVac DOT-Chamber-System!

The layers of the outer cover are welded together dot by dot creating chambers between the dots keeping the internal beads in their correct position, evenly distributed.



- **Optional attachment systems for straps and handles**

Cleats



Hook & Loop




Red Stripes



- **Robust & durable exterior Material**

To achieve outstanding durability, abrasion resistance, tightness, easy cleaning and disinfection, comfort and safety for the patient special materials were chosen. It makes it possible to use the RedVac products in a wide temperature range, from -40°C up to +80°C.



- **Restraint Straps, made of polypropylene (PP) or wipe-clean coated fabric**

Colour coded patient restraint straps, adjustable from both sides, easy to clean, reposition or replace.



- **Vacuum Valve, screw-on, replaceable**

The valve is robust and easy to use, even when wearing gloves. It closes automatically when suction stops preventing loss of vacuum if slipping with the pump. Comes with suitable adapter for most vacuum pumps.



- **Vacuum Valve, welded-on, non-replaceable**

The valve is robust and easy to use, even when wearing gloves. It closes automatically when suction stops preventing loss of vacuum if slipping with the pump.



- **Label System**

An indestructible label with individual serial number and date of manufacture is welded onto every product. Optional labels with property indicator and distributor logo are available.



- **X-ray, CT and MRI compatible**



- **CE**

All RedVac products are CE marked and conform with the European directive 93/42/EEC for medical products class 1.




- **ISO9001 - ISO13485 certified by**




RedVac EMS Vacuum-Mattresses




Body shape

	VM01X02 VM01HX02	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM01CX02 VM01CHX02	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-88-50cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 4		


Body shape

	VM03X03 VM03HX03	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM03CX03 VM03CHX03	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 50-88-50cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 4		



Body shape. Slits for gurney safety belts in neck area.

	VM03X05 VM03HX05	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM03CX05 VM03CHX05	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 50-98-50cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 4		

Body shape. Slits for gurney safety belts in neck area. Glass fiber rod for cervical protection.



	VM03X07 VM03HX07	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM03CX07 VM03CHX07	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 50-88-50cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 4		

Body shape - compact size for helicopter transportation

	VM08X06 VM08HX06	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM08CX06 VM08CHX06	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-78-40cm Length: 195cm Weight: 6kg Carry handles: 8 Patient restraint straps: 4		
	VM09X01 VM09HX01	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM09CX01 VM09CHX01	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-60-40cm Length: 195cm Weight: 5kg Carry handles: 8 Patient restraint straps: 4		

Body shape with head fixation.


Slits for gurney safety belts in neck area. Glass fiber rod for cervical protection.

	VM200X06 VM200HX06	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM200CX06 VM200CHX06	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-88-60cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 6		
	VM201X10 VM201X10	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM201CX10 VM201CHX10	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-100-60cm Length: 210cm Weight: 8kg Carry handles: 8 Patient restraint straps: 6		




Head fixation and slits in neck area for gurney safety belts.

	VM500X03 CLEATS with PP restraint straps VM500HX03 CLEATS with wipe-clean restraint straps VM500CX03 CLICK with PP restraint straps VM500CHX03 CLICK with wipe-clean restraint straps
	Width head/mid/foot: 60-88-88cm Length: 200cm Weight: 7,5kg Carry handles: 8 Patient restraint straps: 5
	VM501X03 CLEATS with PP restraint straps VM501HX03 CLEATS with wipe-clean restraint straps VM501CX03 CLICK with PP restraint straps VM501CHX03 CLICK with wipe-clean restraint straps
	Width head/mid/foot: 60-100-100cm Length: 210cm Weight: 8,5kg Carry handles: 12 Patient restraint straps: 5
	VM502X01 CLEATS with PP restraint straps VM502HX01 CLEATS with wipe-clean restraint straps VM502CX01 CLICK with PP restraint straps VM502CHX01 CLICK with wipe-clean restraint straps
	Width head/mid/foot: 60-80-80cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 5



Head fixation and slits in neck area for gurney safety belts. Glassfiber rod for cervical protection.

	VM5000X01 CLEATS with PP restraint straps VM5000HX01 CLEATS with wipe-clean restraint straps VM5000CX01 CLICK with PP restraint straps VM5000CHX01 CLICK with wipe-clean restraint straps
	Width head/mid/foot: 65-88-88cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 5
	VM5001X01 CLEATS with PP restraint straps VM5001HX01 CLEATS with wipe-clean restraint straps VM5001CX01 CLICK with PP restraint straps VM5001CHX01 CLICK with wipe-clean restraint straps
	Width head/mid/foot: 65-100-100cm Length: 210cm Weight: 8kg Carry handles: 12 Patient restraint straps: 5
	VM5002X01 CLEATS with PP restraint straps VM5002HX01 CLEATS with wipe-clean restraint straps VM5002CX01 CLICK with PP restraint straps VM5002CHX01 CLICK with wipe-clean restraint straps
	Width head/mid/foot: 65-80-80cm Length: 200cm Weight: 6kg Carry handles: 8 Patient restraint straps: 5


Body shape with head fixation. Slits for gurney safety belts. Glass fiber rod for cervical protection. Additional outer protective layer.

	VM7000X03 CLEATS with PP restraint straps VM7000HX03 CLEATS with wipe-clean restraint straps VM7000CX03 CLICK with PP restraint straps VM7000CHX03 CLICK with wipe-clean restraint straps
	Width head-middle-foot: 60-130-60cm Length: 205cm Weight: 10kg Carry handles, integrated in protective layer: 10 Carry handles removable: 4 Patient restraint straps: 5
	VM7001CX04 CLICK with PP restraint straps VM7001CHX04 CLICK with wipe-clean restraint straps
	Width head-middle-foot: 60-110-60cm Length: 205cm Weight: 9,5kg Carry handles, integrated in protective layer: 10 Patient restraint straps: 5
	VM7002CFX01 CLICK with PP restraint straps VM7002CFHX01 CLICK with wipe-clean restraint straps
	Width head-middle-foot: 60-110-60cm Length: 205cm Weight: 9,5kg Carry handles removable: 14 Patient restraint straps: 7




Head fixation. Slits for gurney safety belts in neck area. Glassfiber rod for cervical protection

	<p>VM11X04 CLEATS with PP restraint straps VM11HX04 CLEATS with wipe-clean restraint straps VM11CX04 CLICK with PP restraint straps VM11CHX04 CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-88-80cm Length: 200cm Weight: 7kg Carry handles: 8 Restraint straps: 4</p>
	<p>VM1101X01 CLEATS with PP restraint straps VM1101HX01 CLEATS with wipe-clean restraint straps VM1101CX01 CLICK with PP restraint straps VM1101CHX04 CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-100-80cm Length: 210cm Weight: 8kg Carry handles: 12 Restraint straps: 4</p>
	<p>VM1111X01 CLEATS with PP restraint straps VM1111HX01 CLEATS with wipe-clean restraint straps VM1111CX01 CLICK with PP restraint straps VM1111CHX01 CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-108-80cm Length: 210cm Weight: 8,5kg Carry handles: 12 Restraint straps: 4</p>
	<p>VM1112X01 CLEATS with PP restraint straps VM1112HX01 CLEATS with wipe-clean restraint straps VM1112CX01 CLICK with PP restraint straps VM1112CHX01 CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-108-80cm Length: 220cm Weight: 8,5kg Carry handles: 8 Restraint straps: 4</p>






Rectangular shape. Slits for gurney safety belts in the neck area Glassfiber rod for cervical protection.

	<p>VM11620X02 CLEATS with PP restraint straps VM11620HX02 CLEATS with wipe-clean restraint straps VM11620CX02 CLICK with PP restraint straps VM11620CHX02 CLICK with wipe-clean restraint straps</p> <p>Width: 60cm Length: 200cm Weight: 5kg Carry handles: 8 Restraint straps: 4</p>
	<p>VM11720X01 CLEATS with PP restraint straps VM11720X01 CLEATS with wipe-clean restraint straps VM11720CX01 CLICK with PP restraint straps VM11720CHX01 CLICK with wipe-clean restraint straps</p> <p>Width: 70cm Length: 200cm Weight: 6kg Carry handles: 8 Restraint straps: 4</p>
	<p>VM11820X01 CLEATS with PP restraint straps VM11820X01 CLEATS with wipe-clean restraint straps VM11820CX01 CLICK with PP restraint straps VM11820CHX01 CLICK with wipe-clean restraint straps</p> <p>Width: 80cm Length: 200cm Weight: 7kg Carry handles: 8 Restraint straps: 4</p>
	<p>VM11920X01 CLEATS with PP restraint straps VM11920HX01 CLEATS with wipe-clean restraint straps VM11920CX01 CLICK with PP restraint straps VM11920CHX01 CLICK with wipe-clean restraint straps</p> <p>Width: 90cm Length: 200cm Weight: 7,5kg Carry handles: 8 Restraint straps: 4</p>
	<p>VM1201X01 CLEATS with PP restraint straps VM1201HX01 CLEATS with wipe-clean restraint straps VM1201CX01 CLICK with PP restraint straps VM1201CHX01 CLICK with wipe-clean restraint straps</p> <p>Width: 100cm Length: 210cm Weight: 8,5kg Carry handles: 12 Restraint straps: 4</p>

Tapering shape at head and foot area


	VM10X02 VM10HX02	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM10CX02 VM10CHX02	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 80-88-80cm Length: 200cm Weight: 7kg Carry handles: 8 Restraint straps: 4		
	VM8021X01 VM8021HX01	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM8021CX01 VM8021CHX01	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 80-100-80cm Length: 210cm Weight: 8,5kg Carry handles: 12 Restraint straps: 4		
	VM8121X01 VM8121HX01	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM8121CX01 VM8121CHX01	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 80-108-80cm Length: 220cm Weight: 9kg Carry handles: 12 Restraint straps: 4		

Rectangular shape

	VM620X01 VM620HX01	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM620CX01 VM620CHX01	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
	Width: 60cm Length: 200cm Weight: 6kg Carry handles: 8 Restraint straps: 4	
	VM720X01 VM720HX01	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM720CX01 VM720CHX01	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
	Width: 70cm Length: 200cm Weight: 6,5kg Carry handles: 8 Restraint straps: 4	
	VM820X01 VM820HX01	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM820CX01 VM820CHX01	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width: 80cm Length: 200cm Weight: 7kg Carry handles: 8 Restraint straps: 4		
	VM1021X01 VM1021HX01	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM1021CX01 VM1021CHX01	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width: 100cm Length: 210cm Weight: 8,5kg Carry handles: 12 Restraint straps: 4		
	VM1121X01 VM1121HX01	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	VM1121CX01 VM1121CHX01	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width: 108cm Length: 210cm Weight: 9kg Handles: 12 Carry handles: 12 Restraint straps: 4		

For obese patients


Slits for gurney's safety belts in the neck area. Glass fiber rod for cervical protection

	<p>VM1421X01 CLEATS with PP restraint straps</p>
	<p>VM1421CX01 CLICK with PP restraint straps VM1421CHX01 CLICK with wipe-clean PVC restraint straps</p> <p>Width head-mid-foot: 80-140-140cm Length: 210cm Weight: 11kg Carry handles: 10 Restraint straps: 5</p>

Additional protective layer with integrated handles

	<p>VM2001CX01 CLICK with PP restraint straps VM2001CHX01 CLICK with wipe-clean PVC restraint straps</p> <p>Width head-mid-foot: 80-100-80cm Length: 205cm Weight: 8,5kg Carry handles: 10 Restraint straps: 4</p>
---	---

Additional protective layer with integrated handles



	<p>VM4000CX01 CLICK with PP restraint straps VM4000CHX01 CLICK with wipe-clean PVC restraint straps</p> <p>Width: 95cm Length: 200cm Weight: 8,5kg Carry handles: 10 Restraint straps: 4</p>
---	--

The RedVac Combi-Vacuum-Mattress

2 in 1 - Combines the standard mattress and the vacuum mattress


Body shape with head fixation. Slits for gurney safety belts.

Glass fiber rod for cervical protection.



	<p>Combi-Vacuum-Mattress, without Zipper</p> <p>VM9001CX01 CLICK with PP restraint straps VM9001CHX01 CLICK with wipe-clean PVC restraint straps</p> <p>Width head/middle/foot: 75-104-50cm Length: 200cm Weight: 9kg</p>
	<p>Combi-Vacuum-Mattress, with Zipper</p> <p>VM9002CX02 CLICK with PP restraint straps VM9002CHX02 CLICK with wipe-clean PVC restraint straps</p> <p>Width head/middle/foot: 75-104-50cm Length: 200cm Weight: 9kg</p>
	<p>Vacuum Pillow VK200</p> <p>Width: 57cm Length: 39cm Weight: 0,5kg</p>

The RedVac Vacuum-Head-Immobiliser - for the Spine Board & the Scoop Stretcher

- Quick and easy positioning of the immobiliser
- No manipulation of head and cervical spine
- Adjusts automatically to any head shape and size
- Ears remain free for unrestricted hearing and medical examination
- Solid and comfortable immobilisation
- Vacuum valves: Removable: R-Valve, Fixed: F-Valve



Specification

	<p>Head Immobiliser VK98X07 R-Valve VK98X07V F-Valve Width: 38cm Length: 31cm Weight: 0,7kg</p> <p>Base plates GP08 for scoop stretcher GP11 for back board</p>
	<p>Head Immobiliser with buckles for attachment VK94X17 R-Valve VK94X17V F-Valve Width: 38cm Length: 31cm Weight: 0,7kg</p>

RedVac Vacuum Extrication Device

- the answer to your extrication challenges!
- fits every patient, pregnant – overweight ... no problem!
- comfortable – no respiratory impairment!

Specification

VED24X02	R-Valve
VED24X02V	F-Valve
Width, head: 40cm	
Width, middle: 106cm	
Length: 94cm	
Weight: 2,5kg	



RedVac Vacuum Splints


Optional attachment systems for patient straps:

CLICK HOOK & LOOP RED STRIPES CLEATS



Vacuumvalves:

Removable=R Fixed=F






Patient straps: Washable of PP or Wipe-clean of PVC






RedVac Vacuum Splints "MOZART"







MOZART splints with RED STRIPES & replaceable PP straps

	<p>Forearm Splints (Colour blue – blue)</p> <table border="0"> <tr> <td>VS60X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS60X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS61X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS61X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 3 Width: 35cm Length: 50cm Weight: 0,33kg</p>	VS60X01	Red Stripes, PP-Straps, R-Valve	1-Chamber	VS60X01V	Red Stripes, PP-Straps, F-Valve	1-Chamber	VS61X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	VS61X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber
VS60X01	Red Stripes, PP-Straps, R-Valve	1-Chamber											
VS60X01V	Red Stripes, PP-Straps, F-Valve	1-Chamber											
VS61X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber											
VS61X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber											
	<p>Arm Splints (Colour blue – blue)</p> <table border="0"> <tr> <td>VS62X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS62X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS63X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS63X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width top/foot: 43-70cm Length: 56cm Weight: 0,9kg</p>	VS62X01	Red Stripes, PP-Straps, R-Valve	1-Chamber	VS62X01V	Red Stripes, PP-Straps, F-Valve	1-Chamber	VS63X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	VS63X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber
VS62X01	Red Stripes, PP-Straps, R-Valve	1-Chamber											
VS62X01V	Red Stripes, PP-Straps, F-Valve	1-Chamber											
VS63X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber											
VS63X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber											
	<p>Leg Splints (Colour blue – blue)</p> <table border="0"> <tr> <td>VS64X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS64X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS65X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS65X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width top/foot: 64-70cm Length: 103cm Weight: 1,9kg</p>	VS64X01	Red Stripes, PP-Straps, R-Valve	1-Chamber	VS64X01V	Red Stripes, PP-Straps, F-Valve	1-Chamber	VS65X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	VS65X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber
VS64X01	Red Stripes, PP-Straps, R-Valve	1-Chamber											
VS64X01V	Red Stripes, PP-Straps, F-Valve	1-Chamber											
VS65X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber											
VS65X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber											

MOZART splints with replaceable Loop straps



	<p>Forearm Splints (Colour blue – blue)</p> <table border="0"> <tr> <td>VS60X20V</td> <td>Hook & Loop straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS61X20V</td> <td>Hook & Loop straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 2 Width: 35cm Length: 50cm Weight: 0,33kg</p>	VS60X20V	Hook & Loop straps, F-Valve	1-Chamber	VS61X20V	Hook & Loop straps, F-Valve	Multi-Chamber
VS60X20V	Hook & Loop straps, F-Valve	1-Chamber					
VS61X20V	Hook & Loop straps, F-Valve	Multi-Chamber					
	<p>Arm Splints (Colour blue – blue)</p> <table border="0"> <tr> <td>VS62X20V</td> <td>Hook & Loop straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS63X20V</td> <td>Hook & Loop straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width top/foot: 43-70cm Length: 56cm Weight: 0,9kg</p>	VS62X20V	Hook & Loop straps, F-Valve	1-Chamber	VS63X20V	Hook & Loop straps, F-Valve	Multi-Chamber
VS62X20V	Hook & Loop straps, F-Valve	1-Chamber					
VS63X20V	Hook & Loop straps, F-Valve	Multi-Chamber					
	<p>Leg Splints (Colour blue – blue)</p> <table border="0"> <tr> <td>VS64X20V</td> <td>Hook & Loop straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS65X20V</td> <td>Hook & Loop straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width top/foot: 64-70cm Length: 103cm Weight: 1,9kg</p>	VS64X20V	Hook & Loop straps, F-Valve	1-Chamber	VS65X20V	Hook & Loop straps, F-Valve	Multi-Chamber
VS64X20V	Hook & Loop straps, F-Valve	1-Chamber					
VS65X20V	Hook & Loop straps, F-Valve	Multi-Chamber					

MOZART splints with CLICK & replaceable PP or PVC straps





  	<p>Forearm splints (X03 = red-blue, X05 = blue-blue)</p> <table border="0"> <tr> <td>VS60CX03 / VS60CX05</td> <td>CLICK, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS60CX03V / VS60CX05V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS60CHX03 / VS60CHX05</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS60CHX05V / VS60CHX05V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td colspan="3"> </td> </tr> <tr> <td>VS61CX03 / VS61CX05</td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS61CX03V / VS61CX05V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS61CHX03 / VS61CHX05</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS61CHX05V / VS61CHX05V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 3 Width: 35cm Length: 50cm Weight: 0,33kg</p> <p>Arm splints (X03 = red-blue, X05 = blue-blue)</p> <table border="0"> <tr> <td>VS62CX03 / VS62CX05</td> <td>CLICK, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS62CX03V / VS62CX05V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS62CHX03 / VS62CHX05</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS62CHX05V / VS62CHX05V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td colspan="3"> </td> </tr> <tr> <td>VS63CX03 / VS63CX05</td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS63CX03V / VS63CX05V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS63CHX03 / VS63CHX05</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS63CHX05V / VS63CHX05V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width hand/elbow: 43-70cm Length: 56cm Weight: 0,9kg</p> <p>Leg Splints (X03 = red-blue, X05 = blue-blue)</p> <table border="0"> <tr> <td>VS64CX03 / VS64CX05</td> <td>CLICK, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS64CX03V / VS64CX05V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS64CHX03 / VS64CHX05</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VS64CHX05V / VS64CHX05V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td colspan="3"> </td> </tr> <tr> <td>VS65CX03 / VS65CX05</td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS65CX03V / VS65CX05V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS65CHX03 / VS65CHX05</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VS65CHX05V / VS65CHX05V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width top/foot: 64-70cm Length: 103cm Weight: 1,9kg</p>	VS60CX03 / VS60CX05	CLICK, PP-Straps, R-Valve	1-Chamber	VS60CX03V / VS60CX05V	CLICK, PP-Straps, F-Valve	1-Chamber	VS60CHX03 / VS60CHX05	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber	VS60CHX05V / VS60CHX05V	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber	 			VS61CX03 / VS61CX05	CLICK, PP-Straps, R-Valve	Multi-Chamber	VS61CX03V / VS61CX05V	CLICK, PP-Straps, F-Valve	Multi-Chamber	VS61CHX03 / VS61CHX05	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	VS61CHX05V / VS61CHX05V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber	VS62CX03 / VS62CX05	CLICK, PP-Straps, R-Valve	1-Chamber	VS62CX03V / VS62CX05V	CLICK, PP-Straps, F-Valve	1-Chamber	VS62CHX03 / VS62CHX05	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber	VS62CHX05V / VS62CHX05V	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber	 			VS63CX03 / VS63CX05	CLICK, PP-Straps, R-Valve	Multi-Chamber	VS63CX03V / VS63CX05V	CLICK, PP-Straps, F-Valve	Multi-Chamber	VS63CHX03 / VS63CHX05	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	VS63CHX05V / VS63CHX05V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber	VS64CX03 / VS64CX05	CLICK, PP-Straps, R-Valve	1-Chamber	VS64CX03V / VS64CX05V	CLICK, PP-Straps, F-Valve	1-Chamber	VS64CHX03 / VS64CHX05	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber	VS64CHX05V / VS64CHX05V	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber	 			VS65CX03 / VS65CX05	CLICK, PP-Straps, R-Valve	Multi-Chamber	VS65CX03V / VS65CX05V	CLICK, PP-Straps, F-Valve	Multi-Chamber	VS65CHX03 / VS65CHX05	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	VS65CHX05V / VS65CHX05V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber
VS60CX03 / VS60CX05	CLICK, PP-Straps, R-Valve	1-Chamber																																																																																
VS60CX03V / VS60CX05V	CLICK, PP-Straps, F-Valve	1-Chamber																																																																																
VS60CHX03 / VS60CHX05	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber																																																																																
VS60CHX05V / VS60CHX05V	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber																																																																																
VS61CX03 / VS61CX05	CLICK, PP-Straps, R-Valve	Multi-Chamber																																																																																
VS61CX03V / VS61CX05V	CLICK, PP-Straps, F-Valve	Multi-Chamber																																																																																
VS61CHX03 / VS61CHX05	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																																																																																
VS61CHX05V / VS61CHX05V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																																																																																
VS62CX03 / VS62CX05	CLICK, PP-Straps, R-Valve	1-Chamber																																																																																
VS62CX03V / VS62CX05V	CLICK, PP-Straps, F-Valve	1-Chamber																																																																																
VS62CHX03 / VS62CHX05	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber																																																																																
VS62CHX05V / VS62CHX05V	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber																																																																																
VS63CX03 / VS63CX05	CLICK, PP-Straps, R-Valve	Multi-Chamber																																																																																
VS63CX03V / VS63CX05V	CLICK, PP-Straps, F-Valve	Multi-Chamber																																																																																
VS63CHX03 / VS63CHX05	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																																																																																
VS63CHX05V / VS63CHX05V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																																																																																
VS64CX03 / VS64CX05	CLICK, PP-Straps, R-Valve	1-Chamber																																																																																
VS64CX03V / VS64CX05V	CLICK, PP-Straps, F-Valve	1-Chamber																																																																																
VS64CHX03 / VS64CHX05	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber																																																																																
VS64CHX05V / VS64CHX05V	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber																																																																																
VS65CX03 / VS65CX05	CLICK, PP-Straps, R-Valve	Multi-Chamber																																																																																
VS65CX03V / VS65CX05V	CLICK, PP-Straps, F-Valve	Multi-Chamber																																																																																
VS65CHX03 / VS65CHX05	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																																																																																
VS65CHX05V / VS65CHX05V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																																																																																
  	<p>Forearm splints (X03 = red-blue, X05 = blue-blue)</p> <table border="0"> <tr> <td>VS61CPX03 / VS61CPX05</td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VS61CPX03V / VS61CPX05V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VS61CPHX03 / VS61CPHX05</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VS61CPHX05V / VS61CPHX05V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 3 Width: 35cm Length: 50cm Weight: 0,33kg</p> <p>Arm splints (X03 = red-blue, X05 = blue-blue)</p> <table border="0"> <tr> <td>VS63CPX03 / VS63CPX05</td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VS63CPX03V / VS63CPX05V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VS63CPHX03 / VS63CPHX05</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VS63CPHX05V / VS63CPHX05V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width hand/elbow: 43-70cm Length: 56cm Weight: 0,9kg</p> <p>Leg Splints (X03 = red-blue, X05 = blue-blue)</p> <table border="0"> <tr> <td>VS65CPX03 / VS65CPX05</td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VS65CPX03V / VS65CPX05V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VS65CPHX03 / VS65CPHX05</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VS65CPHX05V / VS65CPHX05V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width top/foot: 64-70cm Length: 103cm Weight: 1,9kg</p>	VS61CPX03 / VS61CPX05	CLICK, PP-Straps, R-Valve	DOT-Chamber	VS61CPX03V / VS61CPX05V	CLICK, PP-Straps, F-Valve	DOT-Chamber	VS61CPHX03 / VS61CPHX05	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	VS61CPHX05V / VS61CPHX05V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber	VS63CPX03 / VS63CPX05	CLICK, PP-Straps, R-Valve	DOT-Chamber	VS63CPX03V / VS63CPX05V	CLICK, PP-Straps, F-Valve	DOT-Chamber	VS63CPHX03 / VS63CPHX05	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	VS63CPHX05V / VS63CPHX05V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber	VS65CPX03 / VS65CPX05	CLICK, PP-Straps, R-Valve	DOT-Chamber	VS65CPX03V / VS65CPX05V	CLICK, PP-Straps, F-Valve	DOT-Chamber	VS65CPHX03 / VS65CPHX05	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	VS65CPHX05V / VS65CPHX05V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																																													
VS61CPX03 / VS61CPX05	CLICK, PP-Straps, R-Valve	DOT-Chamber																																																																																
VS61CPX03V / VS61CPX05V	CLICK, PP-Straps, F-Valve	DOT-Chamber																																																																																
VS61CPHX03 / VS61CPHX05	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																																																																																
VS61CPHX05V / VS61CPHX05V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																																																																																
VS63CPX03 / VS63CPX05	CLICK, PP-Straps, R-Valve	DOT-Chamber																																																																																
VS63CPX03V / VS63CPX05V	CLICK, PP-Straps, F-Valve	DOT-Chamber																																																																																
VS63CPHX03 / VS63CPHX05	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																																																																																
VS63CPHX05V / VS63CPHX05V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																																																																																
VS65CPX03 / VS65CPX05	CLICK, PP-Straps, R-Valve	DOT-Chamber																																																																																
VS65CPX03V / VS65CPX05V	CLICK, PP-Straps, F-Valve	DOT-Chamber																																																																																
VS65CPHX03 / VS65CPHX05	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																																																																																
VS65CPHX05V / VS65CPHX05V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																																																																																

	<p>Arm Splints (Colour blue – blue)</p> <table border="0"> <tbody> <tr> <td>VA14X02</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VA14X02V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VA14HX02</td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VA14HX02V</td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VA14X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14HX01</td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14HX01V</td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </tbody> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p>	VA14X02	Red Stripes, PP-Straps, R-Valve	1-Chamber	VA14X02V	Red Stripes, PP-Straps, F-Valve	1-Chamber	VA14HX02	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber	VA14HX02V	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber	VA14X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	VA14X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber	VA14HX01	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber	VA14HX01V	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber
VA14X02	Red Stripes, PP-Straps, R-Valve	1-Chamber																							
VA14X02V	Red Stripes, PP-Straps, F-Valve	1-Chamber																							
VA14HX02	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber																							
VA14HX02V	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber																							
VA14X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber																							
VA14X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber																							
VA14HX01	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
VA14HX01V	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
	<p>Ankle Splints (Colour blue – blue)</p> <table border="0"> <tbody> <tr> <td>VK15X02</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VK15X02V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VK15HX02</td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VK15HX02V</td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VK15X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VK15X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA15HX01</td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA15HX01V</td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </tbody> </table> <p>Restraint straps: 4 Width: 47cm Length: 57cm Weight: 0,9kg</p>	VK15X02	Red Stripes, PP-Straps, R-Valve	1-Chamber	VK15X02V	Red Stripes, PP-Straps, F-Valve	1-Chamber	VK15HX02	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber	VK15HX02V	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber	VK15X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	VK15X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber	VA15HX01	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber	VA15HX01V	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber
VK15X02	Red Stripes, PP-Straps, R-Valve	1-Chamber																							
VK15X02V	Red Stripes, PP-Straps, F-Valve	1-Chamber																							
VK15HX02	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber																							
VK15HX02V	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber																							
VK15X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber																							
VK15X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber																							
VA15HX01	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
VA15HX01V	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
	<p>Leg Splints</p> <table border="0"> <tbody> <tr> <td>VL21X02</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL21X02V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL21HX02</td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL21HX02V</td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL21X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21HX01</td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21HX01V</td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </tbody> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	VL21X02	Red Stripes, PP-Straps, R-Valve	1-Chamber	VL21X02V	Red Stripes, PP-Straps, F-Valve	1-Chamber	VL21HX02	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber	VL21HX02V	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber	VL21X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	VL21X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber	VL21HX01	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber	VL21HX01V	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber
VL21X02	Red Stripes, PP-Straps, R-Valve	1-Chamber																							
VL21X02V	Red Stripes, PP-Straps, F-Valve	1-Chamber																							
VL21HX02	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber																							
VL21HX02V	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber																							
VL21X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber																							
VL21X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber																							
VL21HX01	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
VL21HX01V	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
	<p>Large Leg Splints (Colour blue – blue)</p> <table border="0"> <tbody> <tr> <td>VL18X02</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL18X02V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL18HX02</td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL18HX02V</td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL18X01</td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL18X01V</td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL18HX01</td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL18HX01V</td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </tbody> </table> <p>Restraint straps: 5 Width: top/hand: 70-45cm Length: 135cm Weight: 2,1kg</p>	VL18X02	Red Stripes, PP-Straps, R-Valve	1-Chamber	VL18X02V	Red Stripes, PP-Straps, F-Valve	1-Chamber	VL18HX02	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber	VL18HX02V	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber	VL18X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	VL18X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber	VL18HX01	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber	VL18HX01V	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber
VL18X02	Red Stripes, PP-Straps, R-Valve	1-Chamber																							
VL18X02V	Red Stripes, PP-Straps, F-Valve	1-Chamber																							
VL18HX02	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber																							
VL18HX02V	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber																							
VL18X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber																							
VL18X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber																							
VL18HX01	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
VL18HX01V	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber																							

VIKING splints with Cleats & replaceable PP or PVC straps

 	<p>Arm Splints (Colour blue – blue)</p> <table border="0"> <tr> <td>VA14X04</td> <td>CLEATS, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VA14X04V</td> <td>CLEATS, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VA14HX04</td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VA14HX04V</td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VA14X03</td> <td>CLEATS, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14X03V</td> <td>CLEATS, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14HX03</td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14HX03V</td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p> <p>Leg Splints (Colour blue – blue)</p> <table border="0"> <tr> <td>VL21X04</td> <td>CLEATS, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL21X04V</td> <td>CLEATS, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL21HX04</td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL21HX04V</td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td>VL21X03</td> <td>CLEATS, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21X03V</td> <td>CLEATS, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21HX31</td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21HX03V</td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	VA14X04	CLEATS, PP-Straps, R-Valve	1-Chamber	VA14X04V	CLEATS, PP-Straps, F-Valve	1-Chamber	VA14HX04	CLEATS, Wipe-clean-Straps, R-Valve	1-Chamber	VA14HX04V	CLEATS, Wipe-clean-Straps, F-Valve	1-Chamber	VA14X03	CLEATS, PP-Straps, R-Valve	Multi-Chamber	VA14X03V	CLEATS, PP-Straps, F-Valve	Multi-Chamber	VA14HX03	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	VA14HX03V	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber	VL21X04	CLEATS, PP-Straps, R-Valve	1-Chamber	VL21X04V	CLEATS, PP-Straps, F-Valve	1-Chamber	VL21HX04	CLEATS, Wipe-clean-Straps, R-Valve	1-Chamber	VL21HX04V	CLEATS, Wipe-clean-Straps, F-Valve	1-Chamber	VL21X03	CLEATS, PP-Straps, R-Valve	Multi-Chamber	VL21X03V	CLEATS, PP-Straps, F-Valve	Multi-Chamber	VL21HX31	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	VL21HX03V	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber
VA14X04	CLEATS, PP-Straps, R-Valve	1-Chamber																																															
VA14X04V	CLEATS, PP-Straps, F-Valve	1-Chamber																																															
VA14HX04	CLEATS, Wipe-clean-Straps, R-Valve	1-Chamber																																															
VA14HX04V	CLEATS, Wipe-clean-Straps, F-Valve	1-Chamber																																															
VA14X03	CLEATS, PP-Straps, R-Valve	Multi-Chamber																																															
VA14X03V	CLEATS, PP-Straps, F-Valve	Multi-Chamber																																															
VA14HX03	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber																																															
VA14HX03V	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																																															
VL21X04	CLEATS, PP-Straps, R-Valve	1-Chamber																																															
VL21X04V	CLEATS, PP-Straps, F-Valve	1-Chamber																																															
VL21HX04	CLEATS, Wipe-clean-Straps, R-Valve	1-Chamber																																															
VL21HX04V	CLEATS, Wipe-clean-Straps, F-Valve	1-Chamber																																															
VL21X03	CLEATS, PP-Straps, R-Valve	Multi-Chamber																																															
VL21X03V	CLEATS, PP-Straps, F-Valve	Multi-Chamber																																															
VL21HX31	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber																																															
VL21HX03V	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																																															
 	<p>Arm Splints (Colour red – blue)</p> <table border="0"> <tr> <td>VA14X07</td> <td>CLEATS, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14X07V</td> <td>CLEATS, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14HX07</td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14HX07V</td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p> <p>Leg Splints (Colour red – blue)</p> <table border="0"> <tr> <td>VL21X07</td> <td>CLEATS, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21X07V</td> <td>CLEATS, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21HX07</td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21HX03V</td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	VA14X07	CLEATS, PP-Straps, R-Valve	Multi-Chamber	VA14X07V	CLEATS, PP-Straps, F-Valve	Multi-Chamber	VA14HX07	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	VA14HX07V	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber	VL21X07	CLEATS, PP-Straps, R-Valve	Multi-Chamber	VL21X07V	CLEATS, PP-Straps, F-Valve	Multi-Chamber	VL21HX07	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	VL21HX03V	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																								
VA14X07	CLEATS, PP-Straps, R-Valve	Multi-Chamber																																															
VA14X07V	CLEATS, PP-Straps, F-Valve	Multi-Chamber																																															
VA14HX07	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber																																															
VA14HX07V	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																																															
VL21X07	CLEATS, PP-Straps, R-Valve	Multi-Chamber																																															
VL21X07V	CLEATS, PP-Straps, F-Valve	Multi-Chamber																																															
VL21HX07	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber																																															
VL21HX03V	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																																															




VIKING splints with CLICK & replaceable PP or PVC straps

 	<p>Arm Splints (X07 = red-blue, X03 = blue-blue)</p> <table border="0"> <tr> <td>VA14CX07 / VA14CX03</td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14CX07V / VA14CX03V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14CHX07 / VA14CHX03</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VA14CHX07V / VA14CHX03V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p> <p>Leg Splints (X07 = red-blue, X03 = blue-blue)</p> <table border="0"> <tr> <td>VL21CX07 / VL21CX03</td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21CX07V / VL21CX03V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21CHX07 / VL21CHX03</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td>VL21CHX07V / VL21CHX03V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	VA14CX07 / VA14CX03	CLICK, PP-Straps, R-Valve	Multi-Chamber	VA14CX07V / VA14CX03V	CLICK, PP-Straps, F-Valve	Multi-Chamber	VA14CHX07 / VA14CHX03	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	VA14CHX07V / VA14CHX03V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber	VL21CX07 / VL21CX03	CLICK, PP-Straps, R-Valve	Multi-Chamber	VL21CX07V / VL21CX03V	CLICK, PP-Straps, F-Valve	Multi-Chamber	VL21CHX07 / VL21CHX03	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	VL21CHX07V / VL21CHX03V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber
VA14CX07 / VA14CX03	CLICK, PP-Straps, R-Valve	Multi-Chamber																							
VA14CX07V / VA14CX03V	CLICK, PP-Straps, F-Valve	Multi-Chamber																							
VA14CHX07 / VA14CHX03	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
VA14CHX07V / VA14CHX03V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
VL21CX07 / VL21CX03	CLICK, PP-Straps, R-Valve	Multi-Chamber																							
VL21CX07V / VL21CX03V	CLICK, PP-Straps, F-Valve	Multi-Chamber																							
VL21CHX07 / VL21CHX03	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
VL21CHX07V / VL21CHX03V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
 	<p>Arm Splints (X07 = red-blue, X03 = blue-blue)</p> <table border="0"> <tr> <td>VA14CPX07 / VA14CPX03</td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VA14CPX07V / VA14CPX03V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VA14CPHX07 / VA14CPHX03</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VA14CPHX07V / VA14CPHX03V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p> <p>Leg Splints</p> <table border="0"> <tr> <td>VL21CPX07 / VL21CPX03</td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VL21CPX07V / VL21CPX03V</td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VL21CPHX07 / VL21CPHX03</td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td>VL21CPHX07V / VL21CPHX03V</td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	VA14CPX07 / VA14CPX03	CLICK, PP-Straps, R-Valve	DOT-Chamber	VA14CPX07V / VA14CPX03V	CLICK, PP-Straps, F-Valve	DOT-Chamber	VA14CPHX07 / VA14CPHX03	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	VA14CPHX07V / VA14CPHX03V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber	VL21CPX07 / VL21CPX03	CLICK, PP-Straps, R-Valve	DOT-Chamber	VL21CPX07V / VL21CPX03V	CLICK, PP-Straps, F-Valve	DOT-Chamber	VL21CPHX07 / VL21CPHX03	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	VL21CPHX07V / VL21CPHX03V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber
VA14CPX07 / VA14CPX03	CLICK, PP-Straps, R-Valve	DOT-Chamber																							
VA14CPX07V / VA14CPX03V	CLICK, PP-Straps, F-Valve	DOT-Chamber																							
VA14CPHX07 / VA14CPHX03	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																							
VA14CPHX07V / VA14CPHX03V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																							
VL21CPX07 / VL21CPX03	CLICK, PP-Straps, R-Valve	DOT-Chamber																							
VL21CPX07V / VL21CPX03V	CLICK, PP-Straps, F-Valve	DOT-Chamber																							
VL21CPHX07 / VL21CPHX03	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																							
VL21CPHX07V / VL21CPHX03V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																							




RedVac Vacuum Splints "GLOBAL"

GLOBAL splints with CLEATS & replaceable Loop straps




7

  	Forearm Splints (Colour blue – blue)		
	VS100X02	CLEATS, Loop straps, R-Valve	1-Chamber
	VS100X02V	CLEATS, Loop straps, F-Valve	1-Chamber
	VS100X01	CLEATS, Loop straps, R-Valve	Multi-Chamber
	VS100X01V	CLEATS, Loop straps, F-Valve	Multi-Chamber
	Restraint straps: 2 Width: 35cm Length: 53cm Weight: 0,5kg		
	Arm Splints (Colour blue – blue)		
	VS120X02	CLEATS, Loop straps, R-Valve	1-Chamber
	VS120X02V	CLEATS, Loop straps, F-Valve	1-Chamber
	VS120X01	CLEATS, Loop straps, R-Valve	Multi-Chamber
	VS120X01V	CLEATS, Loop straps, F-Valve	Multi-Chamber
	Restraint straps: 4 Width: 63-47cm Length: 74cm Weight: 1,1kg		
Leg Splints (Colour blue – blue)			
VS140X02	CLEATS, Loop straps, R-Valve	1-Chamber	
VS140X02V	CLEATS, Loop straps, F-Valve	1-Chamber	
VS140X01	CLEATS, Loop straps, R-Valve	Multi-Chamber	
VS140X01V	CLEATS, Loop straps, F-Valve	Multi-Chamber	
Restraint straps: 6 Width: 70-36cm Length: 110cm Weight: 2,1kg			




GLOBAL splints with Red Stripes & replaceable Loop straps

  	Forearm Splints		
	VS100X04	Red Stripes, Loop straps, R-Valve	1-Chamber
	VS100X04V	Red Stripes, Loop straps, F-Valve	1-Chamber
	VS100X03	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	VS100X03V	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	Restraint straps: 2 Width: 35cm Length: 53cm Weight: 0,5kg		
	Arm Splints		
	VS120X04	Red Stripes, Loop straps, R-Valve	1-Chamber
	VS120X04V	Red Stripes, Loop straps, F-Valve	1-Chamber
	VS120X03	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	VS120X03V	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	Restraint straps: 4 Width: 63-47cm Length: 74cm Weight: 1,1kg		
	Leg Splints		
	VS140X04	Red Stripes, Loop straps, R-Valve	1-Chamber
	VS140X04V	Red Stripes, Loop straps, F-Valve	1-Chamber
	VS140X03	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	VS140X03V	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	Restraint straps: 6 Width: 70-36cm Length: 110cm Weight: 2,1kg		

GLOBAL splints with CLICK & replaceable Loop straps straps



<p>Forearm Splints</p> 	<p>VS100CKX04 VS100CKX04V VS100CKX03 VS100CKX03V VS100CPKX03 VS100CPKX03V</p>	<p>CLICK, Loop straps, R-Valve CLICK, Loop straps, F-Valve CLICK, Loop straps, R-Valve CLICK, Loop straps, F-Valve CLICK, Loop straps, R-Valve CLICK, Loop straps, R-Valve</p>	<p>1-Chamber 1-Chamber Multi-Chamber Multi-Chamber DOT-Chamber DOT-Chamber</p>
<p>Restraint straps: 2 Width: 35cm Length: 53cm Weight: 0,5kg</p>			
<p>Arm Splints</p> 	<p>VS120CKX04 VS120CKX04V VS120CKX03 VS120CKX03V VS120CPKX03 VS120CPKX03V</p>	<p>CLICK, Loop straps, R-Valve CLICK, Loop straps, F-Valve CLICK, Loop straps, R-Valve CLICK, Loop straps, F-Valve CLICK, Loop straps, R-Valve CLICK, Loop straps, R-Valve</p>	<p>1-Chamber 1-Chamber Multi-Chamber Multi-Chamber DOT-Chamber DOT-Chamber</p>
<p>Restraint straps: 4 Width: 63-47cm Length: 74cm Weight: 1,1kg</p>			
<p>Leg Splints</p> 	<p>VS140CKX04 VS140CKX04V VS140CKX03 VS140CKX03V VS140CPKX03 VS140CPKX03V</p>	<p>CLICK, Loop straps, R-Valve CLICK, Loop straps, F-Valve CLICK, Loop straps, R-Valve CLICK, Loop straps, F-Valve CLICK, Loop straps, R-Valve CLICK, Loop straps, R-Valve</p>	<p>1-Chamber 1-Chamber Multi-Chamber Multi-Chamber DOT-Chamber DOT-Chamber</p>
<p>Restraint straps: 6 Width: 70-36cm Length: 110cm Weight: 2,1kg</p>			

GLOBAL splints with CLICK & replaceable PP or PVC straps

<p>Forearm Splints</p> 	<p>VS100CX04 VS100CX04V VS100CHX04 VS100CHX04V VS100CX03 VS100CX03V VS100CHX03 VS100CHX03V VS100CPX03 VS100CPX03V VS100CPHX03 VS100CPHX03V</p>	<p>CLICK, PP-Straps, R-Valve CLICK, PP-Straps, F-Valve CLICK, Wipe-clean-Straps, R-Valve CLICK, Wipe-clean-Straps, F-Valve CLICK, PP-Straps, R-Valve CLICK, PP-Straps, F-Valve CLICK, Wipe-clean-Straps, R-Valve CLICK, Wipe-clean-Straps, F-Valve CLICK, PP-Straps, R-Valve CLICK, PP-Straps, F-Valve CLICK, Wipe-clean-Straps, R-Valve CLICK, Wipe-clean-Straps, F-Valve</p>	<p>1-Chamber 1-Chamber 1-Chamber 1-Chamber Multi-Chamber Multi-Chamber Multi-Chamber Multi-Chamber DOT-Chamber DOT-Chamber DOT-Chamber DOT-Chamber</p>
<p>Restraint straps: 3 Width: 35cm Length: 53cm Weight: 0,5kg</p>			
<p>Arm Splints</p> 	<p>VS120CX04 VS120CX04V VS120CHX04 VS120CHX04V VS120CX03 VS120CX03V VS120CHX03 VS120CHX03V VS120CPX03 VS120CPX03V VS120CPHX03 VS120CPHX03V</p>	<p>CLICK, PP-Straps, R-Valve CLICK, PP-Straps, F-Valve CLICK, Wipe-clean-Straps, R-Valve CLICK, Wipe-clean-Straps, F-Valve CLICK, PP-Straps, R-Valve CLICK, PP-Straps, F-Valve CLICK, Wipe-clean-Straps, R-Valve CLICK, Wipe-clean-Straps, F-Valve CLICK, PP-Straps, R-Valve CLICK, PP-Straps, F-Valve CLICK, Wipe-clean-Straps, R-Valve CLICK, Wipe-clean-Straps, F-Valve</p>	<p>1-Chamber 1-Chamber 1-Chamber 1-Chamber Multi-Chamber Multi-Chamber Multi-Chamber Multi-Chamber DOT-Chamber DOT-Chamber DOT-Chamber DOT-Chamber</p>
<p>Restraint straps: 4 Width: 63-47cm Length: 74cm Weight: 1,1kg</p>			
<p>Leg Splints</p> 	<p>VS140CX04 VS140CX04V VS140CHX04 VS140CHX04V VS140CX03 VS140CX03V VS140CHX03 VS140CHX03V VS140CPX03 VS140CPX03V VS140CPHX03 VS140CPHX03V</p>	<p>CLICK, PP-Straps, R-Valve CLICK, PP-Straps, F-Valve CLICK, Wipe-clean-Straps, R-Valve CLICK, Wipe-clean-Straps, F-Valve CLICK, PP-Straps, R-Valve CLICK, PP-Straps, F-Valve CLICK, Wipe-clean-Straps, R-Valve CLICK, Wipe-clean-Straps, F-Valve CLICK, PP-Straps, R-Valve CLICK, PP-Straps, F-Valve CLICK, Wipe-clean-Straps, R-Valve CLICK, Wipe-clean-Straps, F-Valve</p>	<p>1-Chamber 1-Chamber 1-Chamber 1-Chamber Multi-Chamber Multi-Chamber Multi-Chamber Multi-Chamber DOT-Chamber DOT-Chamber DOT-Chamber DOT-Chamber</p>
<p>Restraint straps: 6 Width: 70-36cm Length: 110cm Weight: 2,1kg</p>			

RedVac Vacuum Splints “REMBRANDT”

REMBRANDT splints with Loop straps



	<p>Arm Splint</p> <p>VS101X02 Hook & Loop straps, R-Valve Multi-Chamber</p> <p>Restraint straps: 2 Width: Top/Hand: 40-30cm Length: 73cm Weight: 0,7kg</p>
	<p>Leg Splint</p> <p>VS103X02 Hook & Loop straps, R-Valve Multi-Chamber</p> <p>Restraint straps: 4 Width: Top/Foot: 60-18cm Length: 96cm Weight: 1,9kg</p>
	<p>Leg Splint</p> <p>VS105X02 Hook & Loop straps, R-Valve Multi-Chamber</p> <p>Restraint straps: 5 Width: Top/Foot: 73-18cm Length: 133cm Weight: 2,2kg</p>

RedVac Vacuum Splints “BASIC”



	<p>Arm Splint</p> <p>VSDS Basic, without straps, F-Valve Width: Top/Hand: 40-30cm Length: 73cm Weight: 0,7kg</p> <p>VSDSX01 Basic with Hook & Loop straps, F-Valve Restraint straps: 3 Width: Top/Hand: 40-30cm Length: 73cm Weight: 0,7kg</p>
	<p>Leg Splint</p> <p>VSDM Basic, without straps, F-Valve Width: Top/Foot: 60-18cm Length: 96cm Weight: 1,9kg</p> <p>VSDMX01 Basic with Hook & Loop straps, F-Valve Restraint straps: 4 Width: Top/Foot: 60-18cm Length: 96cm Weight: 1,9kg</p>
	<p>Leg Splint</p> <p>VSDL Basic, without straps, F-Valve Width: Top/Foot: 73-18cm Length: 133cm Weight: 2,2kg</p> <p>VSDLX01 Basic with Hook & Loop straps, F-Valve Restraint straps: 5 Width: Top/Foot: 73-18cm Length: 133cm Weight: 2,2kg</p>

RedVac Pediatric Vacuum Mattresses

Vacuum Mattresses with CLEATS & PP or PVC straps

	<p> VM710X22 CLEATS, PP-Straps, R-Valve 1-Chamber VM710X22V CLEATS, PP-Straps, F-Valve 1-Chamber VM710HX22 CLEATS, Wipe-clean-Straps, R-Valve 1-Chamber VM710XH22V CLEATS, Wipe-clean-Straps, F-Valv 1-Chamber VM710X08 CLEATS, PP-Straps, R-Valve Multi-Chamber VM710X08V CLEATS, PP-Straps, F-Valve Multi-Chamber VM710HX08 CLEATS, Wipe-clean-Straps, R-Valve Multi-Chamber VM710XH08V CLEATS, Wipe-clean-Straps, F-Valv Multi-Chamber Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 100cm Weight: 3,5kg </p>
	<p> VM710X23 CLEATS, PP-Straps, R-Valve Multi-Chamber VM710X23V CLEATS, PP-Straps, F-Valve Multi-Chamber VM710HX23 CLEATS, Wipe-clean-Straps, R-Valve Multi-Chamber VM710HX23V CLEATS, Wipe-clean-Straps, F-Valv Multi-Chamber Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 100cm Weight: 3,5kg VM712X09 CLEATS, PP-Straps, R-Valve Multi-Chamber VM712X09V CLEATS, PP-Straps, F-Valve Multi-Chamber VM712HX09 CLEATS, Wipe-clean-Straps, R-Valve Multi-Chamber VM712HX09V CLEATS, Wipe-clean-Straps, F-Valv Multi-Chamber Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 100cm Weight: 3,5kg </p>

Vacuum Mattresses with CLICK & PP or PVC straps

 	<p> VM710CX23 CLICK, PP-Straps, R-Valve Multi-Chamber VM710CX23V CLICK, PP-Straps, F-Valve Multi-Chamber VM710CHX23 CLICK, Wipe-clean-Straps, R-Valve Multi-Chamber VM710CHX23V CLICK, Wipe-clean-Straps, F-Valve Multi-Chamber VM710CPX23 CLICK, PP-Straps, R-Valve DOT-Chamber VM710CPX23V CLICK, PP-Straps, F-Valve DOT-Chamber VM710CPHX23 CLICK, Wipe-clean-Straps, R-Valve DOT-Chamber VM710CPHX23V CLICK, Wipe-clean-Straps, F-Valve DOT-Chamber Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 100cm Weight: 3,5kg VM712CX23 CLICK, PP-Straps, R-Valve Multi-Chamber VM712CX23V CLICK, PP-Straps, F-Valve Multi-Chamber VM712CHX23 CLICK, Wipe-clean-Straps, R-Valve Multi-Chamber VM712CHX23V CLICK, Wipe-clean-Straps, F-Valve Multi-Chamber VM712CPX23 CLICK, PP-Straps, R-Valve DOT-Chamber VM712CPX23V CLICK, PP-Straps, F-Valve DOT-Chamber VM712CPHX23 CLICK, Wipe-clean-Straps, R-Valve DOT-Chamber VM712CPHX23V CLICK, Wipe-clean-Straps, F-Valve DOT-Chamber Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 120cm Weight: 3,5kg </p>
---	---

RedVac Pediatric Vacuum Splints

Splints with replaceable Loop straps

Arm splint, angled 	<p> VS56X20V Hook & Loop straps, F-Valve 1-Chamber Restraint straps: 3 Width top/hand: 50-28cm Length: 42cm Weight: 0,6 kg </p>
Leg Splint 	<p> VS58X20V Hook & Loop straps, F-Valve 1-Chamber Restraint straps: 3 Width top/foot: 38-25cm Length: 56cm Weight: 0,7kg </p>

Splints with RED STRIPES & replaceable PP or PVC straps

Arm Splint, angled 	<p> VS56X01 Red Stripes, PP-Straps, R-Valve 1-Chamber VS56X01V Red Stripes, PP-Straps, F-Valve 1-Chamber Restraint straps: 3 Width top/hand: 50-28cm Length: 42cm Weight: 0,6 kg </p>
Leg Splint 	<p> VS58X02 Red Stripes, PP-Straps, R-Valve 1-Chamber VS58X02V Red Stripes, PP-Straps, F-Valve 1-Chamber VS59X01 Red Stripes, PP-Straps, R-Valve Multi-Chamber VS59X01V Red Stripes, PP-Straps, F-Valve Multi-Chamber Restraint straps: 3 Width top/hand: 50-28cm Length: 42cm Weight: 0,6 kg </p>




RedVac Pediatric Restraint System for the ambulance gurney

RedVac JUNIOR

Pediatric Restraint System for use with or without vacuum mattress




JUNIOR is a restraint system for safe transportation of children in ambulances. It is compatible with most ambulance gurneys. The system is suitable for new-borns and children up to around 7 years with body length 130 cm and weight of 35 kg.

<p>Junior</p> 	<p>KF01X07</p> <p>Width head: 39 cm Width foot: 20 cm Length: 74 cm Weight: 1,3 kg</p> <p>Restraint straps: 4 Securing straps: 6</p>		
<p>Vacuum mattresses for JUNIOR</p> 	<p>VKF01PX02 Hook & Loop straps, R-Valve DOT-Chamber VKF01PX02V Hook & Loop straps, F-Valve DOT-Chamber VKF01CPX02 CLICK with PP-Straps, R-Valve DOT-Chamber VKF01CPX02V CLICK with PP-Straps, F-Valve DOT-Chamber VKF01CPHX02 CLICK with Wipe-clean-Straps, R-Valve, R-Valve DOT-Chamber VKF01CPHX02V CLICK with Wipe-clean-Straps, F-Valve DOT-Chamber</p> <p>Width: 45 cm Length: 80 cm Weight: 0,75 kg Restraint straps: 3</p>		
	<p>VKF02PX01 Hook & Loop straps, R-Valve DOT-Chamber VKF02PX01V Hook & Loop straps, F-Valve DOT-Chamber VKF02CPX01 CLICK with PP-Straps, R-Valve DOT-Chamber VKF02CPX01V CLICK with PP-Straps, F-Valve DOT-Chamber VKF02CPHX01 CLICK with Wipe-clean-Straps, R-Valve DOT-Chamber VKF02CPHX01V CLICK with Wipe-clean-Straps, F-Valve DOT-Chamber</p> <p>Width Head/Middle-Foot: 45/62 cm Length: 138 cm Weight: 1,2 kg Restraint straps: 5</p>		

Vacuum Pumps & Carry Bags

	<p>HP01 One-Hand pump, plastic Length: 21cm Width: 14cm Weight: 0,6kg</p>
	<p>HP03 Hand pump, aluminium Length: 30cm Width: 4cm Weight: 0,4kg</p>
	<p>HP04 Hand pump, plastic Length: 28cm Width: 11 / 4,5cm Weight: 0,170kg</p>
	<p>FP01 Foot pump, plastic Length: 43cm Width: 11cm Weight: 0,75kg</p>
	<p>HFP02 Hand- and Foot pump, plastic. Length: 45cm Width: 6cm Weight: 0,95kg</p>
	<p>FP04 Foot pump, aluminium Length: 54cm Width: 6cm Weight: 1,4kg</p>
<p>Carry bag for RedVac Junior</p> 	<p>TS04 Nylon fabric Colour: blue TS04C B/R Cordura fabric Colour: blue, red Wipe-off foil TS04V B/R/G/Y Colour: blue, red, green, yellow Length-Width-Height: 40-20-10cm Weight: 0,25kg</p>
<p>Carry bag for splint set</p> 	<p>TS06 Nylon fabric Colour: blue TS06C B/R Cordura fabric Colour: blue, red Wipe-off foil TS06V B/R/G/Y Colour: blue, red, green, yellow Length-Width-Height: 65-45-12cm Weight: 0,4kg</p>
<p>Carry bag for mattress with backpack straps</p> 	<p>TM01 Nylon fabric Colour: blue TM01C B/R Cordura fabric Colour: blue, red TM01V B/R/G/Y Wipe-off foil Colour: blue, red, green, yellow for mattress with width less than 90cm Length-Width-Height: 80-45-25cm Weight: 0,7kg</p>
	<p>TM02 Nylon fabric Colour: blue TM02C B/R Cordura fabric Colour: blue, red TM02V B/R/G/Y Wipe-off foil Colour: blue, red, green, yellow for mattress with width of more than 90cm Length-Width-Height: 100-45-30cm Weight: 1kg</p>

Accessories for Vacuum Products

<p>Property indicator label</p>		<p>Label sizes: 8cm x 5cm 17cm x 7cm A4</p>
<p>9022 Repair Kit with self-adhesive patches</p>		



TRS Thermo Rescue Sack

for use with a vacuum mattress
and/or spine boards or rescue stretchers

- Water repellent and rugged material provides good protection without additional blanket
- Designed for easy access to the patient
- 10 carry handles
- Under surface with antiskid stripes
- Zipper and Hook-and-Loop closing for fast handling and easy access to the patient
- Optional property indicator label



Undersurface with antiskid stripes

Easy access to the patient

Fits also large Patient

Specification

TRSX01	Material: Nylon, Carry Handles: 10 Width Head/Middle/Foot: 45cm/70cm/40cm Length: 200cm Weight: 2,5kg
TRSX02	Material: Nylon, Carry Handles: 10 Width Head/Middle/Foot: 45cm/70cm/40cm Length: 220cm Weight: 2,9kg
TRSX03	Material: CORDURA quilted, Carry Handles: 10 Width Head/Middle/Foot: 45cm/70cm/40cm Length: 220cm Weight: 2,9kg

HRS Rescue Sack

for use with a vacuum mattress
and/or spine boards or rescue stretchers

- Made of water repellent rugged CORDURA (HRS01) or water repellent heat insulating SOFT SHELL (HRS02)
- Designed for easy access to the patient
- 10 carry handles
- Under surface with antiskid stripes
- Zipper and Hook-and-Loop closing for fast handling and easy access to the patient
- Optional property indicator label
- CE marked



HRS01	Material: Cordura, red, Carry Handles: 10 Width Head/Middle/Foot: 60cm/60cm/40cm Length: 215cm Weight: 2,5kg
HRS02	Material: Softshell, orange, Carry Handles: 10 Width Head/Middle/Foot: 60cm/60cm/40cm Length: 215cm Weight: 2,5kg

BRS Helicopter Rescue Sack

for use with a vacuum mattress
and/or spine boards or rescue stretchers



- Optimised for rescue with helicopter and use with vacuum mattress
- The vacuum mattress can be deflated without opening the rescue bag
- Antiskid under surface to prevent slipping on snow
- Hook and loop fastener and zipper for easy access to the patient and to avoid wind-sack phenomenon and rotation
- Stable head protection
- Lifting strap system is adjustable to achieve horizontal position of the patient
- Patient restraint straps with metal buckles
- Pocket for vacuum pump
- Plated-through holes for free drainage of water
- Anti-rotation rings
- CE labelled and certified for helicopter lifting

BRS03

All straps are sewn onto the rescue bag and are stowed in easy to open pockets



BRS04

The detachable lifting bridle is attached to the carry handles with karabiner



Easy access to the patient



Integrated bag for the vacuum pump



Undersurface with antiskid stripes



Specification



BRS03

BRS04

Carry Handles: 10
Patient restraint straps: 5
Width, Head-Middle-Foot: 40-70-45cm
Length: 220 cm
Weight: 4 kg

TBRS03

Carry bag
Width-Length-Height: 44-80-33 cm
Weight: 1 kg

KOHLBRAT & BUNZ

Emergency & Rescue Systems

Made in Austria

Quality management in accordance with
EN ISO 9001 - EN ISO 13485

certified by



made by

KOHLBRAT & BUNZ GMBH



A-5550 Radstadt

Austria

T: +43 6452 7193 0

F: +43 6452 7193 51

office@kohlbrat-bunz.com

www.kohlbrat-bunz.com

www.redvac.com

Pictures are symbolic and may differ from the original. Stated weights are approximate.
Subject to technical changes.