

TRS – Thermo Rescue Sack **Operating Instructions**



Operating Instructions

Purpose :

The TRS Thermo Rescue Sack is used to protect patients against bad and cold weather during transportation from place of accident to a protected area.

The applications are limited to preclinical assistance and the first phase of clinical support.

Application:

The TRS Rescue Sack is a medical product which may be applied only by medical specialist personnel instructed and trained regularly on the product. The kind of application is in accordance with the current, generally recognised medical guidelines and recommendations of the specialised societies which are concerned.

Cleaning and Care/Disinfection:

The TRS Rescue Sack can be cleaned with luke warm water and a detergent.

Reusability:

The product may be used again if the material is undamaged and its function preserved. Prior to reuse the user must check the proper condition of the product. Only products in perfect technical and hygienic condition may be reused.

Producer :

Kohlbrat & Bunz GmbH, Loretostrasse 4-8, A-5550 Radstadt/Austria

The content of these operating instructions is protected by copyright. Any use, including extracts, shall require written approval.

Subject to technical modifications



Recommended for usage with Vacuum Mattress

- ***10 carry handles***
- ***Undersurface with antiskid stripes***
- ***Hook-and-Loop closing for fast handling and easy access to the patient***
- ***Optional property indicator label***



Specification

TRS01

Color: Red
Carry Handles: 10
Width Head/Middle/Foot: 45cm/70cm/40cm
Length: 200cm
Weight: 2,5kg

TRS02

Color: Red
Carry Handles: 10
Width Head/Middle/Foot: 45cm/70cm/40cm
Length: 220cm
Weight: 2,9kg



**Quality Management System according to
EN ISO 9001:2008 & EN ISO 13485:2003
certified by TÜV Austria**



**KOHLBRAT & BUNZ GMBH
A-5550 Radstadt
Austria
T: +43 (0)6452 7193 0
F: +43 (0)6452 7193 51
office@kohlbrat-bunz.com**