

**EMS Hypothermia  
Prevention Mat**



**Catalog 2024**



**KOHLBRAT  BUNZ**

Emergency & Rescue Systems

Made in Austria

V-01-2024

# Preventing Hypothermia in Patients During Transportation in Cold Environments

- **Robust & durable exterior Material**

To achieve outstanding durability, abrasion resistance, tightness, easy cleaning and disinfection, comfort and safety for the patient special materials were chosen. It makes it possible to use the RedVac products in a wide temperature range, from -40°C up to +80°C.



- **220V Single Charging Port**

For charging one battery at a time



- **220V Dual-Charging Port**

For charging up to two batteries at a time



- **Chargin Cable**



- **Battery Pack**

Can generate temperatures of up to 40C for 2.5 hours



- **Label System**

An indestructible label with individual serial number and date of manufacture is welded onto every product. Optional labels with property indicator and distributor logo are available.



- **X-ray, CT and MRI compatible**



- **CE**

All RedVac products are CE marked and conform with the European directive 93/42/EEC for medical products class 1.



- **ISO9001 - ISO13485 certified by**



# AMW01X01

- *Introducing the AMW01X01 External Heating Unit – an advanced solution designed to elevate patient care during transport in cold environments. Engineered with cutting-edge technology, this heating mat is capable of generating temperatures up to 40°C for approximately 2.5 hours, effectively preventing hypothermia in patients.*
- *The AMW01X01 features a removable battery cleverly housed in a cushioned box at the head area, providing dual functionality by raising the patient's head and ensuring neck stability. Inside the mat, an electrical heating cable dynamically adjusts the temperature based on the surface conditions, offering a reliable and responsive heating solution.*
- *Equipped with three temperature settings, the AMW01X01 adapts to diverse environments and external temperatures, guaranteeing optimal patient comfort and safety. For extended rescue operations, we recommend keeping an extra battery pack on hand, enabling seamless and swift battery replacement to maintain continuous temperature regulation.*
- *With 8 hook loops, versatility is a hallmark of the AMW01X01, allowing seamless integration with various products from our portfolio, including the RollUp stretcher, RedVac mattresses, and the UT2000 stretcher. Its universal compatibility makes the AMW01X01 an indispensable tool for rescue operations, ensuring comprehensive and effective patient care. Elevate your rescue capabilities with the AMW01X01 – where cutting-edge technology meets unparalleled versatility.*



## **Dimonesions:**

Length: 115cm

Width: 43cm

# KOHLBRAT BUNZ

Quality management in accordance with  
EN ISO 9001 - EN ISO 13485

certified by



*made by*

**KOHLBRAT & BUNZ GMBH**

**A-5550 Radstadt**

**Austria**

**T: +43 6452 7193 0**

**F: +43 6452 7193 51**

**office@kohlbrat-bunz.com**

**www.kohlbrat-bunz.com**

**www.redvac.com**

Pictures are symbolic and may differ from the original. Stated weights are approximate.  
Subject to technical changes.